

Cardiology

Defibrillators

Spacelabs Medical
Burdick Products.



DC200 Defibrillator

Specifically, designed for cardiac stress testing or patient-monitored emergency situations. Charges to maximum 360 joules in 9 seconds, Fat and easy 3-step operation, adult or pediatric use, AC line operated with automatic disarm, Weighs less than 17lbs.

BUR161075 DC200



Medic 6 Defibrillator

The compact, easy to use Medic 6 is designed to meet the resuscitation and cardiac therapeutic needs of virtually any setting. Each modality is clearly identified and positioned into dedicated, easy to view areas of the system. Emergency defibrillation buttons are numbered, based on sequence of use,

highlighted and pronounced for fast identification-promoting fast response during critical situations. Lightweight, compact design, 12V rechargeable battery allows for 40+ discharges at 360j or 2+ hours monitoring at full charge. "Precise" Trapezoidal Waveform-using lower peak currents with load compensation.

BUR92450-1A
BUR92450-1AP
BUR92450-1AV

3M Health Care



3M Defib-Pads

Quick and convenient, Safe; eliminates "gel Bridging", years of proven success for defibrillation and cardio versions, 10 pairs/box, 10 boxes/case.

MMM2345N
4-1/2" x 4-1/2"
MMM2345N
4-1/2" x 6"

EKG Machines

Spacelabs Medical
Burdick Products

Eclipse LE II Electrocardiograph

This new version of the original LE model is still the hard-working solution for your basic ECG demands, but now boasts these value-added features:

- One-button automatic or manual operation
- Full-page report-eliminates mounting time and costs
- Intuitive menu functions displayed on an LCD screen
- Full alphanumeric keyboard
- New baseline stabilization filter suppresses baseline wander
- Standard internal rechargeable battery allows the device to be portable and minimizes lines frequency artifact due to poor electrical power
- AHA recommended sampling frequency provides a more accurate analysis-critical for pediatric analysis
- Latest update of the Glasgow ECG analysis program ensures accuracy

BUR92304-1AI09
Eclipse LE II
Interpretive
BUR92304-1AJ09
Eclipse LE II, Non-
Interpretive

BUR92304-1BI09

Eclipse LE II, Interpretive
w/ Stress

BUR92304-1BJ09

Eclipse LE II, Non-
Interpretive, w/ stress



Eclipse 850 & 850i Electrocardiographs

The Eclipse 850 electrocardiograph is the new standard of ECG performance and value, with its compact footprint and large, backlit Three-channel preview display, it is an ideal solution for physician offices, clinics and areas of the hospital where a full-featured system is required but a full-sized system is not needed. Features true, one-button operation, up gradable software, waveform signal quality display, full alphanumeric keyboard automated measurements, optional computerized interpretation, using the Glasgow program, optional storage for 40 ECGs, optional internal modem for fax and SCP-ECG communication

BUR92300-1AL09
Eclipse 850

BUR92300-1BLP09

Eclipse 850 w/ Comm

BUR101368 Eclipse 850

w/Comm/Storage

BUR92300-1AK09

Eclipse 850I w/ Storage

BUR92300-1BK09

Eclipse 850I w/ Comm

BUR101372 Eclipse 850I

w/Comm/Storage

ECG Papers

Hewlett Packard

Z-Fold, Red Grip, Multi Channel

HPM1707A 216mm x

183', 200 sheets/pad, 8

pad/cs

HPM2481A 216mm x

183', 200 sheets/pad, 8

pad/cs

HPM2485A 216mm x

183', 200 sheets/pad, 8

pad/cs

Pad, Red Grip, Multi Channel

HP4700AH-K13

8.5" x 11", 100 sheets/pk,

10pk/cs

HP4700AH-K11

8.5" x 11", 100 sheets/pk,

10pk/cs

Spacelabs Medical

Burdick Products

Single-Channel Mounts

100 mounts/box, 5 box/cs

BUR7857 ECG Card,

pressure sensitive, 12 lead

BUR7909 ECG Card,

pressure sensitive, 3 sec,

12 lead

BUR7037 ECG Card,
pressure sensitive, EK10/
Abridged Elite

BUR7911 ECG Folder,
pressure sensitive, 3 sec,
12 lead

BUR7036 ECG Folder,
pressure sensitive, full
information

ECG Paper

BUR7957 Single

Channel, 45mm x 150ft

roll, black grid, no

header, 10/bx, 100/cs

BUR7940 Single

Channel, 45mm x 150ft

roll, red grid, no header,

10/bx, 100/cs

BUR7958 Single

Channel, 45mm x 150ft

roll, red grid, no header,

10/bx, 100/cs

BUR7975 Single

Channel, 45mm x 150ft

roll, red grid, no header,

10/bx, 100/cs

BUR7932 Single

Channel, 50mm x 150ft

roll, black grid, no

header, 10/bx, 100/cs

BUR7941 Single

Channel, 50mm x 150ft

roll, orange grid, no

header, 10/bx, 100/cs

BUR7006 8.5" x 11" pad,

red grid, no header,

100/box, 10bx/cs

BUR7868 Multi Channel,

8.5" x 11" Z-Fold, red

grid, blank Header,

200/pad, 10pad/cs

BUR7961 Multi Channel,

8.5" x 11", Z-Fold, red

grid, blank header,

200/pad, 10pad/cs

BUR7976 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
200/pad, 10pad/cs
BUR7983 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
200/pad, 10pad/cs
BUR7984 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
200/pad, 10pad/cs
BUR7966 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
200/pad, 10pad/cs
BUR7981 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
300/pad, 10pad/cs
BUR7979 Multi Channel,
8.5" x 11" Z=Fold, red
grid, blank header,
300/pad, 10pad/cs
BUR7989 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
300/pad, 10pad/cs
BUR7988 Multi Channel,
8.5" x 11", Z-Fold, red
grid, blank header,
300/pad, 10pad/cs

ECG Accessories

Spacelabs Medical Burdick Products

ECG Accessories

BUR47036 Adapter Clip,
ECG Universal resting,
Astro trace 10
clips/pouch

BUR7158 Limb strap for
ECG sensor, adult, single
hole, 18"
BUR7518 Patient Cable
ECG Resting, 10 lead
BUR007514/01207000
Patient Cable ECG
Resting 10 Lead
BUR7785 Patient Cable,
ECG Resting 10 Lead

Dynarex Corporation



Electrode Skin Prep Pads

70% isopropyl alcohol
and pumice formulation,
Pumice removes excess
skin oil enhancing
electrodes to skin contact
for EKG testing

DYN1508 100/bx,
10bx/cs

EKG Pads

For fast, easy, convenient
electrode application,
each foil packet contains
five individual cotton
flannel pads premoistened
with a non-irritating
electrolyte solution. Non-
corrosive, will not affect
electrodes. Excellent
conductivity.
DYN1507 5/pk, 50 pk/bx,
10bx/cs

Electrodes

3M Health Care



3M Red Dot Monitoring Electrodes with Foam Tape and Sticky Gel

Foam Tape Backing
Monitoring, can be used
for all cardiac monitoring
applications. The high
adhesion foam backing
makes it ideally suited for
the ER, OR and
diagnostic stress test
applications. High
performance adhesive for
long, secure wear,
includes micro replicated
trace prep surface on
catalog number 2570
version to improve trace
quality. Contains sticky
gel conductor and strip
delivery

MMM2560 Non-Radio
lucent
MMM2570 Radio lucent



3M Red Dot Clear Plastic Monitoring Electrode

3M Red Dot electrodes are designed for hospital wide, clinical situations to minimize the need to remove or replace electrodes when patients move between departments. The clear, plastic back allows for easy viewing of a patient's skin condition throughout the short-term wear cycle. Each electrode has a unique 3M solid gel conductive column and a 3M Sure seal moisture vapor barrier cap to ensure product freshness

MMM2235 Without Abrader, 50/bg, 20bg/cs



3M Red Dot ECG Electrodes

- 3M Sure seal Cap Seal, which prevents dry-out for up to 30 days. Out of the bag
- 3M Solid Gel for non-irritating skin contact that provides clear, sharp and consistent traces
- 3M Built-in Abrader. Clinically proven to lower skin impedance resulting in superior trace quality
- 3M Moistac Adhesive designed for superior, gentle adhesion to any

type of skin with no adhesive residue

Foam Monitoring Electrode, without abrader, 50/bg, 20 bgs/bx, and 1bx/cs.
MMM2237

Foam Monitoring Electrode, with abrader, 25/bg, and 40bg/cs
MMM2259

Monitoring Electrode with Micro pore tape backing, with abrader, 25/bg, 40bg/cs
MMM2249

Soft Cloth Monitoring Electrode, with abrader, 25/bg, 40bg/cs
MMM2255

Wet Gel Monitoring Electrode with Foam backing, with abrader, 25/bg, 40bg/cs
MMM2256

Diaphoretic Foam Electrode, without abrader, 50/bg, 20bgs/bx, 1bx/cs
MMM2330

3M Red Dot Monitoring Electrodes- Alternate Care, 50/bg, 20bgs/bx, 1bx/cs
MMM9640 High Tack Foam tape Backing, pregelled w/ low chloride solid gel

MMM9641 Soft Cloth Ventilated tape Backing w/Hypoallergenic Adhesive, pregelled with Low Chloride Solid Gel

Tyco
Healthcare

Kendall-LTP

Q-Trace 5400 Resting ECG Tab Electrode

The Q-Trace 5400 sets the industry standard for resting ECG tab electrode. This tab electrode features our most popular level of adhesion compatible for most diagnostic applications

LTP31433538 100/bg, 5bg/bx, 4 bx/cs

Q-Trace Gold 5500 Resting ECG Tab Electrode

The Q-Trace Gold 5500 tab electrode features our highest tack adhesive gel and a larger surface area providing maximum adhesion on challenging patients

LTP30807732 100/bg, 5bg/bx, 4bx/cs

Care Resting ECG Tab Electrodes

Available in specialized designs to further address user preferences, CA210

provides a unique, foam substrate design for long term and serial ECG's. CA310 offers rounded corners.

LTPCA210 100/bg,
10bg/bx, 4bx/cs
LTPCA310 100/bg,
10bg/bx, 4bx/cs

Medi-Trace 530 Series Electrodes

Designed for superior performance in Stress, Holter and Diaphoretic applications.

- Conductive Adhesive hydro gel
 - Highly aggressive adhesive
 - Foam substrate conforms easily to skin to ensure electrical contact for consistent tracings
 - Teardrop shape for easy release, removal
 - Pouch packaging contains individual electrodes
- LTP31013926 30/pk,
20pk/cs

IOMED



Trans Q1-GS Electrode

TransQ electrode's unique pigtail connection enables the electrode to

be placed virtually anywhere and to be easily connected to the lead wire. Active Area: 7.6cm², approx. fill Volume: 1.5cc, Maximum Current: 4.0mA, Recommended Dosage: 40mA-minutes, Maximum Dosage: 80mA-minutes, Skin Interface: Gel Sponge, conductive Element: Silver-Silver Chloride

IOM5000017 10
treatment kits per Carton

VERMED



Performance Plus Line

Vermed performance Plus electrodes provide excellent adhesion and fast, clear traces for stress test, holter test and other diagnostic procedures. 5/strip, 5 strip/pouch, 10pouch/bx, 10 bx/cs.

A10005 1-1/2" foam, wet gel
A10006-5 1-7/8" Foam, wet gel
A10006-5T 2-1/4" tape, wet gel



SilveRest Line

SilveRest electrodes are designed for resting EKG and Echo sound procedures. Economical and conveniently

packaged, SilveRest electrodes provide clear traces for accurate diagnosis

A10009 1-3/8" Vinyl tape, wet gel 10/strip, 10 strip/pouch, 10bx/cs
A10013-50S Rect. Foam, wet gel 5/strip, 50 strip/pouch, 20 pouch/cs



Clear Scan line

Clear Scan electrodes are available in various materials, packaging and sizes. These electrodes all have a carbon snap and can be used in cath-lab, MRI, thallium stress and other nuclear medicine procedures. 5/strip, 5 strip/pouch, 100 pouch/cs

A10004-5 1-5/8" Foam, wet gel

Spacelabs Medical Burdick Products



CardioSens ECG Electrodes

The Spacelabs Burdick CardioSens family of resting ECG sensors offer quality you can depend on. Hypoallergenic excellent traces no clean up. Choose from 4 types:

Ultra II, Advantage, K & S.

BUR047029 CardioSens Ultra II
BUR047881 CardioSens Advantage
BUR047866 CardioSens K
BURR047057 CardioSens S



Blue Max Monitoring Electrodes

Blue Max Monitoring Electrodes, used for Holter/Stress/ICU, CCU monitoring, feature a unique offset snap connection which drastically reduces motion artifact while making patient hook-ups a "snap". Blue Max monitoring Electrodes are available in either 12-hour or 50-hour configurations and in a variety of comfortable backings

BUR097001 Vinyl, Snap offset, Wet Gel, 40mm, 25/pouch, 40pouches/cs
BURR097002 Foam, snap offset, Wet Gel, 50mm, 25/pouch, 40pouches/cs



Heart Line Monitoring/Diagnostic Electrodes

Heart Line electrodes provide the highest quality data for virtually any cardiac application. Ideal for long term monitoring and short term diagnostic procedures. 30 electrodes/pouch, 20 pouch/case, case lots only

BUR007784 Soft Cloth Pediatric, 35mm dia., Solid Gel, Monitoring
BURR007778 Comfort tape tear Drop, 37 x 53mm, wet gel, diagnostic
BURR007781 Durafoam Tear drop, 37 x 53mm, wet gel, diagnostic

GELS & CREAM

Parker



Spectra 360 Electrode Gel

The only salt-free and chloride-free electrically conductive gel, recommended for all

electro medical procedures. Salt-free characteristics make it particularly suitable for long-term applications. Spectra 360 differs significantly from all other electrically conductive media...it works by wetting skin resistance

PAR12-08 250g tube
12/bx

Spacelabs Medical Burdick Products



ECG Gel

Cor-Gel

BUR7884 Gel ECG, 150ml tube, 12/bx

Liqui-Cor

BUR7886 Gel ECG, 250ml bottle, 12/box
BUR7887 Gel ECG, 1 liter bottle, 12/box

Weaver & Co



NuPrep Skin Prep Gel

NuPrep gel is intended for use when a reduction of skin impedance would enhance a test result. ECG exams, evoked potential procedures, ECG stress tests, cardiac

rehabilitation monitoring,
cardiac catheter
monitoring exam
procedures are
appropriate for the use of
NuPrep gel
WEA10-30 4oz. tube

Holter Systems

Spacelabs Medical
Burdick Products



Vision Holter Analysis System

Cutting edge technology
for today's cost-conscious
budgets, The Vision
Holter system is a PC-
based system, which
automatically analyzes
and organizes ECG data
from the Holter recorder.
Test results are available
in only minutes.
Meticulous, systematic
artifact identification
makes data review and
editing fast and efficient.
Self-guiding menus, with
on-line videos, make
operation simple even for
the novice. When analysis
is complete, the Vision
Holter system generates a
concise final report as
select by the user

BUR92502-1AHKRQ w/
Cass interface
BUR92502-1AKRQ w/
Recorder
BUR82502-1D w/
Recorder

Spacelabs Medical
Burdick Products

Spacelabs-Burdick Holter Prep Kits

BUR43272 Baseline
Holter, 3 channel
Compact recorder, use
with: 6732 recorder kit
BUR43104 Baseline
Holter, 3 Channel DISC,
use with: 6532 Recorders
Kit
BUR43202 Baseline
Holter, 3 Channel PC
Card Recorder use with:
6632 Recorders Kit
BUR97008 Baseline
Stress Kit

Stress Systems

Spacelabs Medical
Burdick Products.



Extol Stress Testing Systems

Quality and precision in
four easy-to-use
integrated systems, each
system offers

preprogrammed and
custom protocols
including a ramping
protocol, automatic or
manual operation, ST
Segment measurements
and an ultra-quiet, heavy
duty treadmill. Compact
profile fits into limited
space with ease. It is
available in resting
interpretive or non-
interpretive

BUR104265 Extol LE.
BUR104267 Extol Lei
BUR104266 Extol 850
BUR104268 Extol 850i