

## Diagnostic

### PediaPals



#### **PediaPals “Elly” Elephant Otoscope Attachment Diagnostic/Notched Style**


This child-friendly attachment provides a calming effect that makes the visit more pleasant for the child, parent, and healthcare professional. The attachment easily twists on Welch Allyn’s Diagnostic Otoscope. The blue tips (elephant trunk), easily snap on and off.

**PPL100032** 100 per case

#### **PediaPals 2.5mm Blue Otoscope Disposable Tips**

These blue tips (the elephant trunk) easily snap on and off PediaPals “Elly” Elephant Otoscope Attachments. They have a soft, rounded tip for comfort. Overall length is the same as standard tips, with no obstruction of view. Available in bags of 500

**PPL100022** 500/bg,  
10bg/cs



#### **PediaPals “Jamaal” Giraffe Reflex Hammer**

This giraffe easily disarms the child by peeking out of your pocket. Use “Jamaal” to give a child a “Kiss”, instead of a “hit” or “tap”. Ergonomically balanced. Non-toxic hand painted decoration

**PPL100004** 7-1/2

### Welch Allyn



#### **3.5 Halogen Diagnostic Sets**

The portable diagnostic set features Welch Allyn’s finest instruments. Hard cases include molded inserts to keep instruments securely in position, and are made

of attractive, leather-grained black vinyl. The handle has a non-abrasive finish a non-slip grip.

#### **Standard Diagnostic Sets**

**WA97100** Includes  
11710 Ophthalmoscope  
20000 Otoscope, 71000  
handle.

**WA97110** Same as  
WA97100 except with  
26550 Nasal Illuminator

**WA97150** Same as  
WA97100 except with  
25020 Otoscope

#### **Coaxial Diagnostic Sets**

**WA97200** Includes  
11720 Ophthalmoscope  
20000 Otoscope, 71000  
handle.

**WA97200-SM**  
Includes 11720  
Ophthalmoscope  
Diagnostic Otoscope,  
Lithium-Ion (Smart)  
handle

**WA97210** Same as  
WA97200 except with  
26550 Nasal Illuminator

**WA97250** Same as  
WA97200 except with  
25020 Otoscope

#### **Auto Step Diagnostic Set**

**WA97310** Includes  
11730 Ophthalmoscope  
20000 Otoscope, 71000  
handle



**PanOptic  
Ophthalmoscope**

Welch Allyn introduces the PanOptic Ophthalmoscope, a revolutionary new concept for visualizing the fundus. The PanOptic incorporates an innovative and patented new optical design called Axial Point Source Optics. This enabling technology makes it easier to enter small pupils and provides a dramatically wider, more panoramic view of the fundus that's 5X larger than ever before achieved with a standard Ophthalmoscope in an undiluted eye. The ergonomic design of the built-in, soft grip handle provides maximum comfort, balance, and access to controls

WA11810 PanOptic Ophthalmoscope  
WA05258 PanOptic Diagnostic Set Carrying Case



**767 Wall Transformers**

The all-in-one package The 767 Wall Transformer is designed to put everything needed for fast, accurate diagnoses at arm's length. And it does this in a convenient, modular package

WA76710 767 wall Transformer without clock



**Pocket Scope Ophthalmoscope**

Pocket Scope fits into your pocket! Easy to use, lightweight, and perfect for general and specialist use. Halogen light for true tissue color and consistent, long-lasting illumination 12 aperture/filter combinations for greater versatility: micro, small, and large spot sizes,

cobalt blue filter for corneal exams, fixation target, and slit aperture (all on a 6-position aperture dial). Plus: red-free filter and unfiltered halogen on a 2-position filter switch – each position can be engaged with each of the 6 position of the aperture dial. 48 focusing lenses for better resolution, Red-free filter may be used with any aperture.

WA12820 2.5v Pocket Scope Ophthalmoscope with AA Battery Handle  
WA79280 Pocket Scope Desk Charger



**3.5v Halogen Diagnostic Oscope/Throat Illuminator**

Halogen light for true tissue color and consistent, long lasting illumination, Fiber optics for cool light with no reflections, no obstruction. Fibers are flush with distal end for easy cleaning. Built-in throat illumination/penlight simply detaches for high intensity illumination. Wide-angle glass viewing lens combines with large distal opening for

instrumentation under magnification.

WA20000 w/  
Polypropylene Specula,  
sizes 2.5-5mm  
WA03100 Halogen  
replacement Lamp  
WA25020 Same as  
20000 excluding throat  
illuminator



### 3.5v Halogen Operating Otoscope

Open system ideal with  
fully rotatable 5/16"  
diameter magnifying lens  
and rotatable speculum

WA21700 w/5  
Polypropylene Specula  
WA03100 Halogen  
Replacement Lamp



### 3.5v Halogen Pneumatic Otoscope

Halogen light for true  
tissue color and  
consistent, long lasting  
illumination, Rotating  
glass lens for  
instrumentation under  
magnification and ultra-  
seal for superb Otoscope.  
Fiber optic illumination  
for cool light

WA20200 w/ 5  
Polypropylene Specula,  
size 2-9mm  
WA03100 Halogen  
Replacement lamp



### Coaxial Ophthalmoscope

Coaxial optical system  
for ease of entry and  
larger field of view 28  
lens, -25 to +40 diopters  
for excellent resolution.  
Sealed optics, crossed  
linear polarizing filter and  
18 aperture combinations  
for greater versatility.  
Includes Micro spot  
Aperture, rubber brow  
rest, illuminated lens dial,  
and halogen light

WA11720 3.5v  
Coaxial



**Standard  
Ophthalmoscope**  
Six aperture for all  
general and specialist

uses. Micro spot aperture  
allows quick visual entry  
in very small, undilated  
pupils. Rubber brow rest  
prevents scratching of  
eyeglasses. Illuminated  
lens dial clearly identifies  
diopters setting. Halogen  
light for true tissue color  
and consistent, long-  
lasting illumination

WA11710 3.5v  
Standard

### Otoscope Insufflator Bulb and Tube

WA21501 without tip  
WA21501 with Tip



**Kleenspec Disposable  
Specula & Dispenser**  
Single use specula reduce  
risk of cross  
contamination and easily  
attaches to Welch Allyn  
Otoscope. Offers low cost  
per exam. The dispenser  
is a convenient,  
economical dispenser,  
which holds both sizes of  
the 524 series disposable  
tips. Eliminates clutter.

WA52432 Ped. Size,  
1000/bg  
WA52434 Reg. Size,  
1000/bg  
WA52133 3mm,  
500/bg

**WA52134** 4mm,  
500/bg  
**WA52135** 5mm,  
500/bg

**Dispensers**

**WA52100** Dispenser  
(521 series), 10/cs  
**WA52400** Dispenser  
(524 series) 10/cs  
**WA52401** Dispenser  
(524 series), w/ storage  
compartment, 10/cs

**Welch Allyn**



**SofSpec Specula**

Specially contoured for  
ultra pneumatic seal.  
Elastomer tip provides  
more patient comfort.  
May be autoclaved,  
boiled, or sterilized in  
solution. Easily attaches  
to Welch Allyn  
otoscopes.

**For use with all Diagnostic  
Otosopes**

**WA24420** Set,  
includes 3, 5 and 7mm  
Speculum

**Audio Spec Specula**

Set of 3 Audio Spec  
Specula for Audio Scope  
(Small, Medium, Large)  
**WA23400** Set of 3



**Reusable Polypropylene  
Specula**

Great choice for those  
who prefer to use and  
reuse. May be autoclaved,  
boiled, or sterilized in  
solution. Easily attaches  
to Welch Allyn  
otoscopes.

**For use with all Diagnostic  
Otosopes**

**WA24302** 2.5mm  
Speculum  
**WA24303** 3mm  
Speculum  
**WA24304** 4mm  
Speculum  
**WA24305** 5mm  
Speculum  
**WA24400** Set,  
includes 2.5, 3, 4 and  
5mm Speculum

**For use with all  
pneumatic, Operating and  
Consulting Otosopes**

**WA22002** 2mm  
Speculum  
**WA22003** 3mm  
Speculum  
**WA22004** 4mm  
Speculum  
**WA22005** 5mm  
Speculum  
**WA22009** 9mm  
Speculum



**3.5v Rechargeable  
NiCad Handle for  
Desk/Well Chargers**

Lightweight and  
economical desk or well-  
charger handle makes it  
easy and convenient to  
have a fully charged  
handle. Well-balanced  
lightweight and rests  
comfortably in the hand.  
Adjustable rheostat  
controls light intensity.  
Long-lasting rechargeable  
nickel-cadmium, battery  
included. Smooth, non-  
abrasive knurled finish  
ensures a no-slip grip.  
Chrome-plated brass  
construction provides  
durable and reliable long-  
life.


**WA71670** 3.5v  
Rechargeable NiCad  
handle for Desk/Well  
Charges

**3.5v Direct Plug-in  
Rechargeable Handle**

Convenient and  
completely portable, the  
Direct Plug-in will fit all  
Welch Allyn 3.5v  
instrument heads, is  
rechargeable, durable,  
and easy to use.

**WA71000-A** 3.5v direct  
plug-in Rechargeable  
Handle

**WA71050** Bottom  
Section Only



**Welch Allyn  
Replacement Batteries**  
Replacement batteries that will keep your Welch Allyn medical products working for years to come

- WA72200 NiCad Rechargeable Battery (Black)
- WA72300 NiCad Rechargeable Battery (Orange) for 71000-A
- WA72600 2.5v rechargeable battery for Pocket Scope
- WA72700 7.2v Rechargeable battery for Original Micro Tympanometer
- WA72000 2.5v Nickel-Cadmium Rechargeable battery (Red)
- WA72100 2.5v Rechargeable Battery (green)
- WA72500 7.2v Rechargeable Battery (2)

## Anoscopes

### Cooper Surgical



**Sani-Scope Anoscope**

Clear plastic material allows 360° visibility and reduces thermal shock. Sani-Scope is disposable which allows protection for the doctor, the nurse and the patient from cross-contamination. Single use obturator simplifies insertion and allows little discomfort to patient. 25 units/bgs, 4 bgs/cs

CPS82420 Ones Size, Anoscope

## Welch Allyn

### Kleenspec Disposable Anoscope with Obturator

A high-performance disposable, anoscope for convenient safe and efficient examination with the finest illumination available, Individually wrapped for safety and convenience. 25/bxs, 4 bx/cs

WA53110 Kleenspec Disposable Anoscope  
WA73305 Transformer with 5-ft. Cord  
WA36019 Fiber Optic Light head

## Lighting

### AARON

Medical Industries, Inc.



### Fluoro-Dot Penlight

This AA-size light for use with fluorescein dye is excellent for diagnosing abrasions and fitting contact lenses. Features include, a special written No.47A cobalt blue filter, two-year shelf life and four-hour burn time.

AA6802 Fluoro-Dot

### Fluoro-Slit Penlight

Designed with a special bulb and Wratten 47A cobalt blue filter to project a concentrated cobalt-slit image for diagnostic use with fluorescent dye. The Fluoro-Slit has a two-year shelf life and a four-hour continuous burn time

AA6803 Fluoro-Slit

### Cobalt Blue Penlight

Disposable penlights have a 1 year shelf life and are available with a cobalt blue filter for ophthalmic procedures. Also available in replaceable battery and lamp format

AA6303 Penlight, 3/pk

## BURTON



### **Burton Ultraviolet (Woods) Lights**

For Physician's offices, emergency rooms, and outpatient facilities. Burton UV Lights are universally accepted for diagnosing disease and injuries to the skin and eye, including detection of corneal scratches, foreign bodies in the eye, blocked tear ducts, fungus infections, Erythasma, Tinea Vesicolor, Vitiligo, Scabies and Microsporium Canis. Portable hand held lights for maximum flexibility. Optically ground and polished magnifier lens with 3 power magnification and focal length of 8"

BUT0316030 Magnifier and 2 white bulbs and 2 UV bulbs, 115v

## **Welch Allyn**



### **Direct Focusing Headlight**

Lightweight and comfortable to wear, Convenient focusing sleeve for adjusting the

spot size according to the procedure. "Quick Detach" electrical connector provides mobility from room to room

WA46003 Direct Focusing Headlight Complete



**Halogen Headlight**  
Halogen illumination in a lightweight and portable headlight for the office exam room, Provides true tissue color and consistent, long-lasting illumination. Iris diaphragm regulates diameter of light spot. "Quick Detach" connector on headband for greater mobility

WA49003 Halogen Headlight Complete  
WA02600 6v Halogen Lamp

## **Vaginal Speculas**

## **Cooper Surgical**



### **Sani-Spec Disposable Specula**

Sani-Spec Disposable Specula provide all the benefits of stainless steel specula without the risk of cross-contamination. Cost efficient and easy to use. Allows the examination of vaginal mucosa at 360 degree, angle without rotation. Adjusts for elevation and angle and holds in place. Made from hard, clear, non-toxic plastic. 10/bgs, 10bg/cs

CPS82329 Small  
CPS82330 Medium  
CPS82331 Large

**Sani-Spec Spec Light**  
Gyne-light with adapter, reusable with AAA batteries

CPS82355 Gyne-Light with adapter

## **Welch Allyn**



**Kleenspec 580 Series  
Disposable Vaginal  
Specula**

Halogen light for true tissue color and consistent, long lasting illumination Disposable for economy, safety and convenience One-handed operation with angulation and elevation for convenience especially narrow blade width available for adolescent and geriatric patients. Individually wrapped. Graves style. 4/cs

- WA58000** Small
- WA58001** Medium
- WA58004** Large
- WA78103** Complete Permanent Illumination System
- WA78010** Illumination System w/ switch
- WA78600** Expendable Illumination for Vaginal Specula, 6/bx, 6bx/cs

**KleenSpec 586 Series  
Disposable Vaginal  
Specula**

Halogen light for true tissue color and consistent, long lasting illumination Disposable

for economy, safety and convenience One-hand operation with angulation and elevation for convenience Expendable light carrier for safety. Pederson style.  
WA58600 Small, 25/bx, 4bx/cs  
WA58601 Medium, 25/bx, 4bx/cs

**X-Ray Supplies**

**Image Plus**

**Image Plus X-Ray Film**

**Rare Blue Earth**

- DID10533 14x17
- DID10690 30x35
- DID10532 14x36
- DID10536 11x14
- DID10537 14x36
- DID10689 7x17
- DID10534 10x12
- DID10538 24x30
- DID10539 24x24
- DID10535 8x10
- DID10540 18x24

**Blue**

- DID10549 14x17
- DID10553 30x35
- DID10543 14x36
- DID10548 11x14
- DID10550 14x36
- DID10545 7x17
- DID10547 10x12
- DID10552 24x30
- DID10551 24x24
- DID10546 8x10
- DID10544 18x24

**Duplicating**

- DID10557 14x17
- DID10555 10x12

**Green**

- DID10521 14x17
- DID10529 30x35

- DID10691 14x36
- DID10525 11x14
- DID10526 14x36
- DID10523 7x17
- DID10524 10x12
- DID10522 24x30
- DID10528 24x24
- DID10520 8x10
- DID10527 18x24

**Graham-Field**

**X-Ray Illuminators**

Designed for desk viewing or wall hanging. Each utilizes 15W fluorescent tubes for even diffusion of light. Housings are stainless steel with baked enamel finish. Features instant grip film retainers and Plexiglas viewing panel. All units come with three-wire grounded cords.

**Two Bank**

Two 14" x 17" viewing areas in a single case, each bank has own on/off switch

- GF3793GF 1ea

**Kimberly-Clark  
Corporation**

**X-Ray Gown**

- KC69766 Large, Blue, 100/cs

**BANTA  
Healthcare Group**



### **X-Ray Storage/Mailing Envelopes**

X-Ray Envelopes may be used for record storage or hospital filing

**TID950221** Storage, Envelope, Kraft, 10-1/2" x 12-1/2", 100/cs

**TID950218** Storage, Heavyweight, Kraft, 14-1/2" x 17-1/2", 500/cs

**TID950220** Storage, Lightweight, Kraft, 14-1/2" x 17-1/2", 500/cs

**TID950219** Mailer Envelope w/ insert, Kraft, 15" x 18", 50/cs