

Defibrillators

MINIDEF 3

The easy, reliable and efficient defibrillator

The handling of the MINIDEF 3 defibrillator is very safe and easy. After switching it on, the device performs an internal self-test with error-detection automatically. The clearly structured operating and display elements inform the user unmistakably about the operating state of the defibrillator.

DEFIGARD 1002

The compact defibrillator with unique monitoring capabilities

The DEFIGARD 1002 is a defibrillator with monitor, recorder and unique modular monitoring capabilities for pre-hospital and hospital use. Inserting the modules into the defibrillator provides instant on-screen results without prior configuration. Interchangeable adapters allow the use of external or internal electrodes, adhesive pads or optional "fibrillation recognition" in automated external defibrillation (AED) mode.

DEFIGARD 3002

The in-hospital defibrillator with all defibrillation modes

A set of 4 patient cable cartridges provides for:

- External defibrillation with hand-held paddles
- External defibrillation with large single use self-adhesive electrodes
- Internal defibrillation with "spoon" electrodes
- Semi Automated External Defibrillation (SAED) with optional VF/VT detection

DEFIGARD 6002

DEFIGARD 6002 - the top of the range defibrillator

Features:

- ECG
- SpO₂
- NIBP
- FICO₂
- Defibrillator
- Temperature
- Battery and mains connection
- Color screen
- External Pacemaker
- Alarms
- Trends
- 3-channel recorder
- Mains and battery operated, incl. charger
- Recorder providing 12-lead ECG and SpO₂ print-outs with associated measurements
- Memory for 500 events and trends over the past 2 or 6 hours (ECG, SpO₂)



Defibrillators


THE DEFIBRILLATION LINE

Emergencies do not tolerate any compromises


Modular defibrillator equipment. Conceived for medical emergencies.
Considering our product philosophy: "Keep it simple and save".

The convincing advantages of FRED®, FRED® easy and FRED® easyport:


- Wide operating range (hospitals, established doctors, rescue units, public access)
- Pictograms and self-explanatory program flow for secure handling
- Automatic functions prevent mal-operations
- Miscellaneous options such as ECG, SpO₂, electronic and wireless data transfer etc.



FRED®



FRED® easy



The innovation of the year 2003:

FRED® easyport!

- Light: only 450 grams, incl. batteries
- Small: only 133 x 126 x 35 mm
- High-resolution LCD

Multipulse Biowave®

SCHILLER takes a significant leap in cardiac defibrillation. The biphasic technology of **Multipulse Biowave®** opens new perspectives for the rate of survival after out-of-hospital cardiac arrest. With SCHILLER's dedication to save lives, a new generation of transthoracic defibrillation has started!

This new Defibrillation technology is being used for the first time worldwide in the semi-automatic **FRED®** (First Responder External Defibrillator), **FRED® easy** and **FRED® easyport** from SCHILLER.

