

Resting ECG

CARDIOVIT AT-101

High-tech ECG in a compact, budget-priced device


Performance:

- Intuitive, easy use with direct function keys
- 12-channel Resting ECG
- Graphical LC display
- Storage and serial data transfer for Resting ECG data on PC
- Various formatted printouts on internal thermal printer, external DeskJet or - Laser printer
- Memory: Saving possibility for 40 records (optional)
- Measurements
- Interpretation (optional)
- Thrombolysis (optional)

Features:

- Integrated 3-channel monitor
- Alphanumeric keyboard
- Direct keys for the most important functions
- Easy Print: Integrated high-resolution thermal printer, subsequent printout on thermoreactive paper possible
- Interface for external printer for A4 printouts on regular paper
- Interface for data-transmission

Possibility of data transfer to SCHILLER Data Management System SEMA-200.



MINISCOPE MS-3

Emergency 1-channel ECG with built-in tripod electrode

Features:

- Display time of 46 seconds
- Memory for up to 30 minutes of ECG recording
- 3 or 5 lead patient cable connectable
- PC program for evaluation/automatic analysis
- Printout through any ECG unit
- Weighs only 250 g

Pocket ECG which can be used at any time & anywhere!

Also available for veterinary medicine.

